



Sciex 6500+ QTrap Operational Steps for Trained Personnel

Last Updated 09172017


1. If any of the following instructions does not make sense to you or was not covered during your hand-on training, stop immediately and inform IMSERC staff
2. If Analyst is running, double left click on the Hardware Configuration, deactivate the hardware profile, then exit the program, and click on yes to all save changes requests



3. Run Services (double left click on the icon on the desktop)
 - a. Highlight AnalystDeviceDriveService and then click on **Stop** the service
 - b. Now click on AnalystService and then click on **Stop** the service



Stop the service
Restart the service



4. Install your column (ensure you have capped both ends of the one you remove)
5. Check mobile phase volumes for all the lines to ensure that they are above 40% fill level and that you have enough for your entire run set
6. You may have to coordinate with IMSERC staff to refill the bottles in the beginning or the middle of your run
7. Double Click on Analyst icon on the desk top





8. Select your project in Analyst
 - a. Your project will be set up for you during your hands-on training.




9. Activate the Hardware Profile that corresponds to your method



- a. The green arrow indicates the profile was activated properly
- b. The divert valve on the MS goes to A (source) position, hold down the button under the LED light for a split second to switch it back to B (waste)



10. Click on AD icon on the bottom menu bar (not the one on the Analyst menu) to bring up ADD window
 - a. Click on Status





11. You should see Multisampler in green and Binary Pump and Column Compartment in brown color banners
12. Hover the cursor over each box and click on the small resize icon to maximize the modules information panels




13. Do the following for each module
 - a. Multisampler

- i. Auto-Clean; this primes the needle wash port. You should see wash solvent bubble out of the wash port
- ii. Prime Solvents; check to see that solvent drips out of each line, a and b
- iii. Assign Wellplates to all locations that you will be using (only approved plates that have a valid plate map programmed into the software by IMSERC staff can be used; if you cannot pick a map from the menu, then you are not allowed to use the plate)
- iv. Factory default maps have a star at each end of their names
- v. Those programmed by IMSERC staff have no stars
- vi. If you are not completely certain that you are using the right plate name, stop and ask IMSERC staff for help.
- vii. Wrong plate name can lead to damaging the inject needle assembly and you will be responsible for repair cost and downtime.




b. Binary Pump

- i. Check and if needed switch solvent selection valve A and B
- ii. Prepare Pump
- iii. Seal Wash Prime




c. Column Compartment

- i. You can change the temperature set point if it is significantly different that the temperature requested in your method.
- ii. You can change it manually by a right click on the method item and entering the new value



14. If you have not already, ensure that you are in your project folder



15. Create a new or open an existing acquisition batch

- a. Enter a title for the SET.
 - i. This will be the name of your *.wiff file that will keep the data generated with every injection. Individual injections do not get their own data file name as is typical of other platforms, but instead, they are kept in a *.wiff file.
 - ii. Click on Add Set.
 - iii. A batch file can have more than one set. Each set could be a different variable that you are studying. For example, we could have a set called new column and another called old column, in a hypothetical scenario where were comparing two columns with the same set of samples
- b. Select the method you want to use from the drop down menu
 - i. During your training, a method will be written for your project
 - ii. You will learn how to make minor changes to your method and save it with your name attached
 - iii. The steps required for running a method will be part of a more extensive training
 - iv. IMSERC staff will work with your on an ongoing basis to deliver the advanced training and it may require more than one session
- c. Click on Add Samples and if you don't know how many samples, just enter 12 to get started

The screenshot shows the 'Select Method for Sample Set' dialog box with the following settings:

- Set: New Column Test - 09132017
- Quantitation: none
- Acquisition: Caffeine-LC/MSM-08212017
- Use as Template:
- Use Multiple Methods:

Below the dialog box is a table with the following columns: Sample Name, Rack Code, Rack Position, Plate Code, Plate Position, Vial Position, Data File, and Inj. Volume (µl).

Sample Name	Rack Code	Rack Position	Plate Code	Plate Position	Vial Position	Data File	Inj. Volume (µl)
1 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
2 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
3 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
4 Caffeine-1ng-mL	Multi Drawer	1	54ViaGenPlateC	3	2	New Column Test - 09132017	1.000
5 Caffeine-10ng-mL	Multi Drawer	1	54ViaGenPlateC	3	3	New Column Test - 09132017	1.000
6 Caffeine-25ng-mL	Multi Drawer	1	54ViaGenPlateC	3	4	New Column Test - 09132017	1.000
7 Caffeine-100ng-mL	Multi Drawer	1	54ViaGenPlateC	3	5	New Column Test - 09132017	1.000
8 Caffeine-250ng-mL	Multi Drawer	1	54ViaGenPlateC	3	6	New Column Test - 09132017	1.000
9 Caffeine-500ng-mL	Multi Drawer	1	54ViaGenPlateC	3	7	New Column Test - 09132017	1.000
10 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
11 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
12 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
13 Sample001	Multi Drawer	1	54ViaGenPlateC	3	8	New Column Test - 09132017	1.000
14 Sample002	Multi Drawer	1	54ViaGenPlateC	3	9	New Column Test - 09132017	1.000
15 Sample003	Multi Drawer	1	54ViaGenPlateC	3	10	New Column Test - 09132017	1.000
16 Sample004	Multi Drawer	1	54ViaGenPlateC	3	11	New Column Test - 09132017	1.000
17 Sample005	Multi Drawer	1	54ViaGenPlateC	3	12	New Column Test - 09132017	1.000
18 Sample006	Multi Drawer	1	54ViaGenPlateC	3	13	New Column Test - 09132017	1.000
19 Sample007	Multi Drawer	1	54ViaGenPlateC	3	14	New Column Test - 09132017	1.000
20 Sample008	Multi Drawer	1	54ViaGenPlateC	3	15	New Column Test - 09132017	1.000
21 Sample009	Multi Drawer	1	54ViaGenPlateC	3	16	New Column Test - 09132017	1.000
22 Sample010	Multi Drawer	1	54ViaGenPlateC	3	17	New Column Test - 09132017	1.000
23 Sample011	Multi Drawer	1	54ViaGenPlateC	3	18	New Column Test - 09132017	1.000
24 Sample012	Multi Drawer	1	54ViaGenPlateC	3	19	New Column Test - 09132017	1.000
25 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
26 QC Std 1	Multi Drawer	1	54ViaGenPlateC	3	20	New Column Test - 09132017	1.000
27 Blank	Multi Drawer	1	54ViaGenPlateC	3	1	New Column Test - 09132017	1.000
28 Sample013	Multi Drawer	1	54ViaGenPlateC	3	21	New Column Test - 09132017	1.000
29 Sample014	Multi Drawer	1	54ViaGenPlateC	3	22	New Column Test - 09132017	1.000
30 Sample015	Multi Drawer	1	54ViaGenPlateC	3	23	New Column Test - 09132017	1.000
31 Sample016	Multi Drawer	1	54ViaGenPlateC	3	24	New Column Test - 09132017	1.000

- i. Sample Name
 1. Fill the sample name column
 2. You must run three blanks to get the UPLC column conditioned
 3. Fill the rest of the column with the sequence of samples you want to run
 4. A typical run list is blanks, calibration curve, blanks, 12 samples, blank, known standard, blank, 12 samples, etc.
- ii. Rack Code
 1. MultiDrawer
- iii. Rack Position
 1. Always 1
- iv. Plate code
 1. Click in the cell and from the drop down menu, select *54VialGenPlateC for standard 2 mL LC vials
 2. Click in the cell and from the drop down menu, select *96VialGenPlateC for 96 wellplate (shallow or deep)
 3. Note that specific plate maps are assigned when preparing the Multisampler (section 13a.iii)
- v. Plate Position
 1. From the drop down menu in each cell, select number 1 through 8 (since we have installed full height drawers) where the samples have been loaded
- vi. Vial Position

1. Numbers will be 1 through 54 or 1 through 96, depending on the actual plates being used
 2. You will work with multiples of 6 for 54 vial plate and multiples of 8 for the 96 well wellplates
 3. For 54 vial plates, vials 1 through 6 are located at A1, B1, through F1, and then number seven is A2 location, etc.
 4. For 96 well wellplate, positions 1 through 8 are at A1, B1, through H1, and then number 9 is at A2 location
 5. After you have loaded the samples, you install the plate by rotating it 90 degrees counter clockwise so that the A1 position is now at front-left
- vii. Data File
1. This column should be populated with the Set name
 2. You can override it by simply typing a different title in each cell
- viii. Inj. Volume (uL)
1. This value is imported automatically from the method you selected.
 2. You can override it by typing a different number
 3. You must make manual note of the volume change since the data file information, after the data is collected, will only show you the volume that was entered in the selected method and not the override value

d. Click on Locations tab

- i. Click on the tray where the samples are loaded, example shown is position 3.
- ii. Plate position numbering starts at the bottom front for 1, bottom back for 2, one level higher and front is 3, and so on ...
- iii. At the top left corner, click on the forth icon, the square filled with red dots
- iv. You will now see the location of each vial with a map of the plate



1. Red circles will have the acquisition row number in the middle
2. Green circles means you are injecting from that vial location more than once
3. Hover the cursor on one of the circles and vial number will be revealed

16. Once you finished filling out batch file sample tab, save it by clicking on File on the top menu, and then click on save as.

- a. Your naming convention should include Your Last Name – PI Last Name – Date
- b. You can have this name applied to other files as long as file extensions are different

17. If the compounds that are being analyzed have a quantitation method already (see below), then you can assign a quantitation method and fill out the quantitation tab

18. Click on submit tab



- If any errors are found, they will show in the submit status box
- For example, if you use multiple Rack Codes, you will have to correct it in the Sample tab before you can submit samples
- Once all errors are addressed, you can highlight each line in any order and submit, i.e. to be added to the queue
- Analyst allows you to continuously submit samples
- If you submit samples that you may not want to run, go to the queue tab and highlight and delete the submission(s)

Sample Name	Rack Position	Plate Position	Vial Position	Acquisition Method	Quantitation	Data File	Set Name	Submit Status
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Caffeine-10ng-m	1	3	2	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Caffeine-10ng-l	1	3	3	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Caffeine-50ng-l	1	3	4	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Caffeine-100ng-l	1	3	5	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Caffeine-250ng-l	1	3	6	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Caffeine-500ng-l	1	3	7	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample001	1	3	8	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample002	1	3	9	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample003	1	3	10	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample004	1	3	11	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample005	1	3	12	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample006	1	3	13	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample007	1	3	14	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample008	1	3	15	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample009	1	3	16	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample010	1	3	17	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample011	1	3	18	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample012	1	3	19	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
QC-std-1	1	3	20	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Blank	1	3	1	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample013	1	3	21	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample014	1	3	22	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample015	1	3	23	Caffeine-LCMRM-0821	none	New Column T	New Column	Not
Sample016	1	3	24	Caffeine-LCMRM-0821	none	New Column T	New Column	Not

19. Click on queue

- View queue icon is on the top left
- All of the submitted samples should be in a table with the status Waiting
- Ready & Equilibrate
 - Click on the icon with a red check mark, to the right of the horizontal hourglass, to get the instrument ready; all three hardware modules should be green
 - Click on the icon without a check mark, to the left of the horizontal hourglass, to equilibrate the column.
 - If you had turned on all of the modules through the ADD, and the mobile phase for the initial condition had been running for more than 5 minutes, you can skip the equilibrate step

Start Time	Sample Name	Plate Pos	Vial Pos	Status
9/13/2017 1:35:23 PM	Blank	5	54	Acquired
9/13/2017 1:40:21 PM	Blank	5	54	Acquired
9/13/2017 1:44:59 PM	Blank	5	54	Acquired
9/13/2017 1:49:38 PM	Blank	5	54	Acquired



20. Click on Start Sample

- Hover the mouse over the first Erlenmeyer flask from the left
- When all conditions are satisfied and the system is ready to inject the first sample, the Start Sample Erlenmeyer flask will have a solid black color, otherwise it will be gray

21. Viewing Data

- Click on Explore
- Double Click on Open Data File
- Select Sample window opens
- Select the *.wiff file in the Data Files box that corresponds to your data
- Click on a sample in the Samples box
- The data corresponding to that injection is now displayed
- If you want to navigate to other data, use the black arrows
- Right clicking on the experimental title will open up the scan types, the assigned names and the transition table, in the case of MRM experiments

Select Sample dialog box:



Please select a data file and a corresponding sample:

Data Files	Samples
20170509-US-prot-phosph.wiff	20170509-US-prot-phosph-1
29-12-2016	20170509-US-prot-phosph-1000
Acetaminophen_20170306181211	20170509-US-prot-phosph-100
Acetaminophen_20170306181359	
Ag-Cell-Assay-1.wiff	
B-Gal-nano-06272017.wiff	
Bupropion-cmpd-opt_2017022621015	
Caffeine.wiff	
Christina-08092017.wiff	
Christina-08102017.wiff	
Christina-E1-08102017.wiff	
Christina-E-strogens-08092017.wiff	


Analyst software interface showing the 'Explore Mode' and various data visualization options.

22. Quantitate


- a. Click on the quantitate segment
- b. Build Quantitation Method
 - i. Double click on Build Quantitation Method
 - ii. In the Select Sample window, select the *.wiff file that contains the representative sample data
 - iii. Select a representative sample, one that you are sure the detector is not saturated or the S/N is not acceptable – for example a standard with a concentration in the middle of the linear dynamic range






- iv. Click on Okay
- v. Now in the Components tab, assign the transitions to the Internal Standards and the Analytes tables
- vi. You can only use values available in the drop down menus, as you click on each cell
- vii. Don't forget to remove the internal standard transition from the Analytes table and fill in the Internal Standard column in the Analytes table



- viii. Click on the Integration tab
- ix. Inspect the peak detected
- x. You can change the peak finding criteria (RT window and expected RT, for example) if the peak selected is not the desired peak
- xi. To include peak tailing/fronting, change the peak width value
- xii. You can also change the units of the concentration values
- xiii. Click on Calibration tab



- xiv. Typically, the fit is linear and no weighting is use and the conditions applied to all analytes
 - xv. While inspecting the results table, any changes to these conditions and/or excluding data points from the standard curve can be tested
 - xvi. Now save the quantitation method by clicking on file and perform a save as operation
 - xvii. The saved file will have a *.qmf extension
- c. Quantitation Wizard
- i. Double click on quantitation wizard
 - ii. Select a *.wiff file and highlight all of the samples that you want to include in the quantitation processing and click on the right arrow to transfer them to the Selected Samples box
 - iii. Then click on next twice
 - iv. Select the method you wrote for this set from the drop down menu
 - v. Click Finish
 - vi. A results table will now be shown
 - vii. Save this table
 - viii. Click on the calibration pane to see how good the calibration fit is
 - ix. During the training you will learn how to exclude points and discuss acceptance criteria
 - x. You can also open and view other saved results table by clicking on Review Results Table



23. Reporting Data

- Click on Companion Software
- Double click on the Reporter 3.2
- Try different report formats and save the report as a word or pdf file

