


# T1 Experiment

## Introduction

When an NMR sample sits in the magnet, the applied static magnetic field  $B_0$  will generate the equilibrium magnetization  $M_0$  along +z axis. When a RF pulse is allied to the sample, the net magnetization will be rotated away from +z axis. T1 relaxation (longitudinal or spin-lattice) is the process by which the net magnetization goes back to its initial maximum value ( $M_{z,eq}$ ) parallel to  $B_0$ .

The inversion-recovery experiment measures  $T_1$  relaxation times of any nucleus. If the net magnetization is placed along the -z axis, it will gradually return to its equilibrium position along the +z axis at a rate governed by  $T_1$ . The equation governing this behavior as a function of the time  $t$  after its displacement is:

$$M_z(t) = M_{z,eq} \left( 1 - 2e^{-t/T_1} \right)$$




The basic pulse sequence consists of an  $180^\circ$  pulse that inverts the magnetization to the -z axis. During the following delay, relaxation along the longitudinal plane takes place. Magnetization comes back to the original equilibrium z-magnetization. A  $90^\circ$  pulse creates transverse magnetization. The experiment is repeated for a series of delay values taken from a variable delay list. A 1D spectrum is obtained for each value of  $vd$  and stored in a pseudo 2D dataset. The longer the recycle delay ( $d1$ ) is, the

more precise the  $T_1$  measurement is. Ideally  $d1$  should be set to  $5 \cdot T_1$ . A rough estimation of the  $T_1$  value can be calculated from the null-point value by using  $T_1 = t_{null} / \ln 2$ .

## Setting up proton T1 Experiment

- 1) To set up a  $T_1$  experiment, start with recording a normal proton spectrum to adjust the spectral sweep width **SWH**, acquisition time **aq** and other parameter if necessary.
- 2) Create new dataset and load "**Proton\_T1**" parameter set. Update the parameters with the ones you obtained from last step. The recycle delay **D1** should be  $\sim 2 \cdot 5 \cdot T_1$ . Adjust **NS** accordingly to give sufficient S/N (**fig 1**).



- 3) Edit the "t1delay" by clicking on  at VDLIST line in **fig 1**. **Fig 2** is a good starting list.
- 4) Change the "TD" value for F1 dimension to the number in your VDLIST (**fig 3**)
- 5) Collect the pseudo 2D  $T_1$  dataset

**Fig 1. ACQUPARS display in “pulse program parameters” view**

The screenshot shows the ACQUPARS window with the 'General' tab selected. The title bar indicates 'Probe: CP QCI 600S3 H/F-C/N-D-05 Z'. The 'General Channel f1' section contains the following parameters:

Parameter	Value	Description
PULPROG	t1lr	Pulse program for acquisition
TD	32786	Time domain size
SWH [Hz, ppm]	8196.72	Sweep width
AQ [sec]	1.9999460	Acquisition time
RG	64	Receiver gain
DW [µsec]	61.000	Dwell time
DE [µsec]	20.00	Pre-scan-delay
D1 [sec]	5.000000000	Relaxation delay; 1-5 * T1
d11 [sec]	0.0299999993	Delay for disk I/O [30 msec]
DS	0	Number of dummy scans
NS	2	Scans to execute
VDLIST	t1delay	Variable delay list
vd [sec]	5.00000000	vd[10]={ 5.000000 sec 0.001000 sec... }

**Fig 2. An example of t1delay list with 8 delays**

The screenshot shows a text editor window titled 't1delay (/opt/topspin4.0.8/exp/stan/nmr/lists/vd)'. The list contains the following values:


Line	Value
1	0.050
2	0.100
3	0.250
4	0.500
5	0.800
6	1.5
7	3
8	5


**Fig 3. ACQUPARS display in “all acquisition parameters” view**

The screenshot shows the ACQUPARS window with the 'Experiment' tab selected. The title bar indicates 'Probe: CP QCI 600S3 H/F-C/N-D-05 Z'. The 'Experiment' section contains the following parameters:


Parameter	Value	Description
PULPROG	t1lr	Current pulse program
AQ_mod	DQD	Acquisition mode
FnTYPE	traditional(planes)	nD acquisition mode for 3D etc.
FnMODE	QF	Acquisition mode for 2D, 3D etc.
TD	32786	Size of fid
DS	0	Number of dummy scans
NS	2	Number of scans
TD0	1	Loop count for 'td0'
TDav	0	Average loop counter for nD experiments

## Processing


- 1) Process and adjust phase for the dataset. Use **rser n** (n is the number of total delays) to read out the last fid. Process and phase correct it. On the Adjust Phase toolbar, click **Save for spectrum**. 
- 2) Go back to pseudo 2D T1 dataset by closing the 1D window
- 3) At the command prompt, type **xf2** to process only the F2 axis. Type **abs2** to baseline correct the rows.
- 4) On the menu bar, click **Applications**.
- 5) On the **Dynamics** button, click the drop-down arrow to see more options and in the list, select **T1/T2 Module**.




- 6) The flow buttons change to determine the T1 / T2 relaxation times. While executing the steps below, message windows will be displayed. Please read each message thoroughly and follow the instructions. On the Workflow button bar, click **Fid**



- 7) In the Extract a row from 2d data window, click **Spectrum**




- 8) Enter Slice Number = n (the last one).




- 9) On the Workflow button bar, click **Peaks/Ranges**  Peaks/Ranges

- 10) In the Define Peaks and/or Integrals window, click **Manual Integration**.




11) Define the regions by drawing an integral over the peaks of interest, On the Integration toolbar, click **Save/export integration regions** 

12) In the list, select **Export Region To Relaxation Module**.



13) In the Prepare relaxation data window, click **OK**



14) On the Workflow button bar, select **Relaxation**.




15) By default, the selected areas are peak-picked, and the first peak is displayed in the Relaxation window.

16) On the Workflow button bar, select **Fitting**





17) In the message window, click **Close**.



18) In the Relaxation parameters window, click **OK** and select **Area** as Fitting type.

19) On the Workflow button bar, select **Calculation**. > Calculation

20) In the message window, click **Close**.



21) In the T1/T2 tools bar, click **Calculate fit for all peaks** >>



22) On the Workflow button bar, select **Report** Report

