

FIVE steps to use IMSERC NMR

1. Login to [NUCORE](#) with your **netid** and logon to the instrument:
 - a) turn on the computer monitor at instrument,
 - b) your usage count starts
2. Login to instrument with your **operator id** (usually same as netid)
3. Load your sample and run your experiment
4. Logout from the instrument
5. Login to [NUCORE](#) and logout your instrument session:
 - a) turn off the computer monitor at instrument,
 - b) your usage count stops

Common Commands/Parameters in TopSpin


Setting up experiments & Processing

sx 10	to put #10 (could be any position) sample on the autosampler into magnet.
sx ej	to put sample inside magnet back to the autosampler
rga	automatically set receiver gain
zg	start acquisition
tr	transfer data (while acquisition is in progress)
multizg	start multiple acquisitions starting from current dataset
go	submit experiment to acquisition
stop	abort an acquisition, losing all the FID data recorded so far
halt	halt the running acquisition, saving the recorded FID data to hard disk
efp	weighted Fourier Transformation for 1D dataset
apk	do automatic phase correction
abs	automatically optimize baseline
xfb	weighted Fourier Transformation for 2D dataset

Important parameters for acquisition

P1	F1 channel 90° pulse width, micro seconds
P2	F1 channel 180° pulse width
RG	Receiver gain
D1	relaxation delay, 1 to 5 times T1
2TD	Time domain data points for F2 (direct dimension)
2SW	spectral width in ppm for F2 (direct dimension)
1TD	Time domain data points for F1 (indirect dimension)
1SW	spectral width in ppm for F1 (indirect dimension)
AQ	Acquisition time in seconds
NS	Number of scans
DS	Number of dummy scans
NUC1 – NUC8	Nucleus observed (1H, 13C, 31P, 19F, etc.)
O1 – O8	Frequency offset for channel 1 – 8 in Hz
O1P – O8P	Freq. offset for channels 1 – 8 in ppm
SFO1 – SFO8	Freq. for channels 1 – 8 in MHz

Workflow to interactively setup experiment with TopSpin



1. Create new dataset and setup initial parameters

2. Insert/Eject Sample

3. Lock on your selected Solvent

4. Tune/Match the Probe to the nucleus of your expt

5. Shim

6. Load probe related parameters

7. Auto set receiver gain

8. Acquire fid

9. Process data

The red arrowed fields below need to be filled, red circled parameters in next page can be changed based on your need

- Sample name
- Experiment number
- Experiment to run
- Solvent used
- Directory for the dataset (your folder under your group name)

The screenshot shows the 'Create New Dataset - new' dialog box. The fields are filled with the following values:

- NAME: yw1-042
- EXPNO: 3
- PROCNO: 1
- Experiment: 13C
- Options: Set solvent (C2D2Cl4), Execute 'getprosol'
- DIR: /home/walkon/data/Marks/ygq763
- Number of additional datasets: 1

Parameters setting for 1D and 2D experiments

1 presat 1 1 /home/walkon/data/zhang

Spectrum ProcPars **AcquPars** Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Probe: PA BBO 600S3 BB-H-D-05 Z BTO

General Channel f1

General

PULPROG	zgpr	...	E	Pulse program for acquisition
TD	32768			Time domain size
SWH [Hz, ppm]	9615.38	16.0212		Sweep width
AQ [sec]	1.7039360			Acquisition time
RG	203			Receiver gain
DW [µsec]	52.000			Dwell time
DE [µsec]	6.50			Pre-scan-delay
D1 [sec]	2.000000000			Relaxation delay; 1-5 * T1
d12 [sec]	0.00002000			Delay for power switching [20 usec]
DS	0			Number of dummy scans
NS	4			1 * n, total number of scans: NS * TD0
TD0	1			Number of averages in 1D

Channel f1

SFO1 [MHz]	600.1678208			Frequency of ch. 1
O1 [Hz, ppm]	2820.78	4.700		Frequency of ch. 1
NUC1	1H	Edit...		Nucleus for channel 1
P1 [µsec]	12.700			F1 channel - 90 degree high power pulse
PLW1 [W, dB]	31.623	-15.00		F1 channel - power level for pulse (default)
PLW9 [W, dB]	0.00020402	36.90		F1 channel - power level for presaturation

1 tp35 12 1 /home/walkon/data/zhang

Spectrum ProcPars **AcquPars** Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Probe: PA BBO 600S3 BB-H-D-05 Z BTO

Experiment Width Receiver Nucleus Durations Power Program Probe Lists NUS Wobble Lock Automation Miscellaneous User Routing

Experiment

PULPROG	hsqcetgp	...	E	Current pulse program
AQ_mod	DQD			Acquisition mode
FnTYPE	traditional(planes)			nD acquisition mode for 3D etc.
FnMODE		Echo-Antiecho		Acquisition mode for 2D, 3D etc.
TD	1024	256		Size of fid
DS	16			Number of dummy scans
NS	8			Number of scans
TD0	1			Loop count for 'td0'
TDav	0			Average loop counter for nD experiments

Width

SW [ppm]	20.0264	165.0000		Spectral width
SWH [Hz]	12019.230	24902.283		Spectral width

2 RR 3 1 /home/walkon/data/zhang

Spectrum ProcPars **AcquPars** Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Probe: PA BBO 600S3 BB-H-D-05 Z BTO




Experiment Width Receiver Nucleus Durations Power Program

Nucleus 1

NUC1	1H	Edit...	13C	Observe nucleus
O1 [Hz]	2820.78		11318.37	Transmitter frequency offset
O1P [ppm]	4.700		75.000	Transmitter frequency offset
SFO1 [MHz]	600.1678208		150.9229277	Transmitter frequency

Setup experiment using ICONNMR w/ SampleXpress

1. Login with your operator ID
2. Load your sample to SampleExpress
3. Click the Holder # where you sample is loaded.
4. Click Add
5. Fill in following fields: Name, No., Solvent, and Experiment
6. Change parameters if needed
7. Click Submit
8. Logout ICON by clicking on “change user”



Setup experiment using ICONNMR w/o autosampler

1. Login with your operator ID
2. Click the Holder 1
3. Click Add
4. Fill in following fields: Name, No., Solvent, and Experiment
5. Change parameters if desired (add additional experiments if needed)
6. Click Submit
7. Insert your sample when promoted
8. (optional) If you sample is not ejected automatically after finishing, switch lock sample in

